

QUAT SANITIZER



- EPA Approved
- Tabletops
- Hard Surfaces
- Counters
- Prep Areas
- Cook Areas
- Equipment
- Free Test Strips
- No Down Time
Waiting for
Tablets to
Dissolve
- 180 Pillow Packs
Equals 6 Jars of
Tablets



PACK SIZE:

180 / 9 gram
quick dissolving
pillow packs

*Store in a dry place.

MADE IN USA

ChemPro™

PROFESSIONAL

1-800-227-7765

- Quat Sanitizer is an EPA Approved Third Sink Sanitizer that can also be used as a table top sanitizer.
- Pillow packs *dissolve quickly* and are pre-measured which helps ensure that the sanitizer is at proper strength and that the strict guidelines and directions that pertain to this EPA Approved Product are met.
- Properly stored, Quat Sanitizer can maintain a strength of *200ppm for 2 - 3 months.*

*Follow instructions to make 200ppm Quat Sanitizer.
Quat test strips included.*

USE AS A THIRD SINK SANITIZER:

MIXING INSTRUCTIONS:

- 1) Fill sink with water.
- 2) Add 1 quick dissolving pillow pack of Quat Sanitizer per every 5 gallons of water.
- 3) Allow product to mix well.
- 4) Soak items with sanitizer and allow to remain wet for at least 60 seconds.
- 5) Rack to dry. Do not rinse.

USE TO SANITIZE ITEMS; SUCH AS, TABLES, COUNTERS, AND COOKING AREAS:

MIXING INSTRUCTIONS:

- 1) Place 1 (one) quick dissolving pillow pack into the mixing tote labeled Quat Sanitizer.
- 2) Add 5 gallons of water.
- 3) Allow product to mix well.

USE INSTRUCTIONS:

- 1) Fill directly from tote valve into your labeled, color coded, spray bottle.
- 2) Spray area to soak with sanitizer and allow to remain wet for at least 60 seconds.
- 3) Drain and air dry. Do not rinse.

300 N. Lake St. • Montgomery, IL 60538

ultimate
cleaning
system